


WIDEX HEARING AIDS - PERFORMANCE & FEATURE OVERVIEW

FEATURES	BENEFITS	MOMENT 440	MOMENT 330	MOMENT 220	MOMENT 110	MAGNIFY 110	MAGNIFY 50
PureSound™ Program	Provides an easily recognised and transparent sound for first time users with mild hearing loss (RIC and BTE devices only)	✓	✓	✓	✓		
TruAcoustics™	Provides a more natural and personalised sound for each individual ear	✓	✓	✓	✓	✓	✓
SoundSense Learn	An App feature which lets you select your preferred sound experience anywhere, anytime.	mRIC RD, RIC312D, BTE13D, BTE RD	mRIC RD, RIC312D, BTE13D, BTE RD	mRIC RD, RIC312D, BTE13D, BTE RD			
SoundSense Adapt	Trains the feature sets to your listening needs via an onboard control, DEX or App	✓	✓	✓			
Wind Noise Attenuation	Better speech understanding and comfort in windy/ outdoor conditions.	✓					
High Frequency Boost	Improved music appreciation, and ability to hear high pitched sounds e.g. birds.	✓					
Impact Program	Personalised program for clients who want a sharper, crisper overall sound.	✓					
Speech Enhancer	Best speech understanding in noise through accurate and instant enhancement of speech to focus on the dominant voice in the crowd.	Real Time InterEar	InterEar				
Digital Pinna	Provides natural sound in quiet environments by directing sounds into the ear.	✓	✓				
Trusound Softener	Ensures sudden sounds are heard comfortably, and are not jarring or harsh.	✓	✓	✓			
Directional Focus	Enables you to focus on speech from behind (onboard program) or from behind or either side (App function). Useful in situations like the car.	✓	✓	✓	✓	✓	✓
Made for iPhone & Android	Streams phone calls and audio directly from a smartphone to your hearing aid.	mRIC RD, RIC312D, BTE13D, BTE RD	mRIC RD, RIC312D, BTE13D, BTE RD	mRIC RD, RIC312D, BTE13D, BTE RD	mRIC RD, RIC312D, BTE13D, BTE RD	mRIC RD, RIC312D, BTE13D, BTE RD	RIC312D, BTE13D
Noise Reduction	Reduces background noise to enhance listening comfort.	✓	✓	✓	✓	✓	✓
Soft Level Noise Reduction	Improves listening comfort by reducing sounds like a humming fridge or fan.	✓	✓	✓	✓	✓	✓
True Input Technology	Improved speech understanding and sound quality, without distortion, even in very loud situations.	✓	✓	✓	✓	✓	✓
Feedback Cancelling	Uninterrupted listening - no annoying whistling from the hearing aid.	InterEar	InterEar	InterEar	✓	✓	✓
Zen+ (InterEar)	Multiple programs of Soothing chimes and tones which help you manage your tinnitus.	✓	✓	✓	✓	✓	✓
SmartSpeak	Aid will tell you when you need to change battery or tell you the program you are on.	✓	✓	✓	✓	✓	✓
Audibility Extender	A solution for dead regions of the cochlear, to hear important speech sounds like “S”.	✓	✓	✓	✓	✓	
Sound classes	Automatically detects your listening environment.	11 InterEar	7 InterEar	4	3	3	1
Processing channels	For a more customised fitting to your audiogram.	15 InterEar	12 InterEar	10	6	6	4
High Definition Locator	Finer and better isolation of speech when it's noisy to focus on the speech in front of you and reduce the background noise.	15	12	10	Broadband	Broadband	Broadband



SOUND CLASSES



QUIET



QUIET
with speech



MUSIC
classical



MUSIC
modern



PARTY



PARTY
with speech



TRANSPORT



TRANSPORT
with speech



URBAN



URBAN
with speech



SOCIAL



UNIVERSAL



PURE SOUND



ZEN

AUTOMATIC PROGRAMS

TINNITUS

										
MODEL	mRIC R D	RIC 312 D	RIC 10	BTE R D	BTE 312	BTE 13 D	ITE/ITC	XP	CIC	CIC-MICRO
STYLE	RIC / RITE			BTE		POWER BTE	ITE / ITC		CIC	
RECEIVERS / TUBES	V2M, V2P	V1 S,M,P,HP	V1 S,M	V2 Thin Tubes or Earhook	V1 Thin Tubes or Earhook	V2 Thin Tubes or Earhook				
MOMENT 440/330/220/110	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MAGNIFY 110	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
MAGNIFY 50		✓	✓		✓	✓		✓	✓	✓
Compatible with CROS-Fusion RIC and CROS-Fashion BTE transmitters			✓		✓		✓	✓	✓	
POWER SOURCE										
Rechargeable	Li-Ion			Li-Ion						
Battery Size		312	10		312	13	312	312	10	10
Battery Life hrs <i>(varies with streaming)</i>	16-20	155-165	100	24-37	155	230	125-135	135-150	105	105
Fitting Range dB <i>(varies with receiver)</i>	0-110	0-115	0-90	0-105	0-95	0-105	0-90	0-85	0-90	0-90
USER CONTROLS										
Separate VC/ program				✓		✓	option			
Configurable button	✓	✓		yes + LED indicator	✓					
RC-DEX control	✓	✓	✓	✓	✓	✓	✓	✓	✓	
App control	MOMENT™ or MAGNIFY™	MOMENT™ or MAGNIFY™	TONELINK™	MOMENT™ or MAGNIFY™	TONELINK™	MOMENT™ or MAGNIFY™	TONELINK™	TONELINK™	TONELINK™	TONELINK™
Telecoil		✓		✓		✓		✓		
AUDIO STREAMING										
2.4 audio streaming (iOS)	✓	✓		✓		✓				
Android audio streaming (ASHA)	✓	✓		✓		✓				
TV-PLAY	✓	✓		✓		✓				
TV-DEX, COM-DEX, COM-DEX Remote Mic	✓	✓	✓	✓	✓	✓	✓	✓	✓	

EASYWEAR™ INSTANT EAR-TIPS

V1 EASYWEAR S, M, P RECEIVERS & V2 EASYWEAR M & P RECEIVERS

V1 & V2 EASYWEAR THIN TUBES (0.9 & 1.4 MM)

OPEN SIZE: XS, S, L

ROUND SIZE: S, M, L

DOUBLE SIZE: XS, S, M, L

TULIP SIZE: M, L









ONE-VENT/TWO-VENT

CUSTOM EAR-TIPS

V1 EASYWEAR S, M, P RECEIVERS & V2 M, P RECEIVERS

V1 & V2 EASYWEAR THIN TUBES (0.9 & 1.4 MM)


HARD (HOLLOW)


HARD EXTENDED (SOLID)

SOFT EXTENDED (SOLID)

HARD EMBEDDED (ENCASED)









TRANSPARENT/BEIGE
Variable venting options

Ability to fit smaller ears with larger receivers.

No ability to change receiver/ear wire in clinic. Full remake required.

Want a more in-depth explanation of Widex features and products? Visit Widex Learn at <http://learn.widex.pro/learn> for interactive learning modules.