

EARPIECE OVERVIEW

V.2 RIC MODELS

HEARING LOSS	MODELS/FAMILY	RECEIVER	RIC/RITE	CUSTOM/INSTANT
MILD TO MODERATE	mRIC R D WIDEX MOMENT™	V.2 M	V.2 EASYWEAR WIRED RIC Size: 0, 1, 2, 3, 4	INSTANT OPEN EAR-TIP XS, S, L TULIP EAR-TIP M, L DOUBLE EAR-TIP XS, S, M, L ROUND EAR-TIP, one-vent/two-vent; S, M, L CUSTOM OPEN EAR-TIP, TRANSPARENT/BEIGE HARD EAR-TIP, TRANSPARENT/BEIGE SOFT EAR-TIP, TRANSPARENT HARD/SOFT EXTENDED EAR-TIP, TRANSPARENT
SEVERE	mRIC R D WIDEX MOMENT™	V.2 P	V.2 EASYWEAR WIRED RIC Size: 0, 1, 2, 3, 4	INSTANT TULIP EAR-TIP M, L DOUBLE EAR-TIP XS, S, M, L ROUND EAR-TIP, one-vent/two-vent; S, M, L CUSTOM HARD EAR-TIP, TRANSPARENT/BEIGE SOFT EAR-TIP, TRANSPARENT HARD/SOFT EXTENDED EAR-TIP, TRANSPARENT

V.2 THIN TUBES

HEARING LOSS	MODELS/FAMILY	TUBE	RIC/RITE	CUSTOM/INSTANT
MILD TO SEVERE	BTE 13 D WIDEX MOMENT™ , EVOKE™	V.2 THIN TUBE 0.9 V.2 THIN TUBE 1.4	V.2 EASYWEAR THIN TUBE Size: 0, 1, 2, 3, 4	INSTANT OPEN EAR-TIP XS, S, L TULIP EAR-TIP M, L DOUBLE EAR-TIP XS, S, M, L ROUND EAR-TIP, one-vent/two-vent; S, M, L CUSTOM OPEN EAR-TIP, TRANSPARENT/BEIGE HARD EAR-TIP, TRANSPARENT/BEIGE SOFT EAR-TIP, TRANSPARENT HARD/SOFT EXTENDED EAR-TIP, TRANSPARENT

EARPIECE OVERVIEW

RIC/RITE MODELS

HEARING LOSS	MODELS/FAMILY	RECEIVER	RIC/RITE	CUSTOM/INSTANT
MILD	RIC 10 WIDEX MOMENT™ , EVOKE™ FUSION 2 WIDEX EVOKE™ FUSION WIDEX EVOKE™ , UNIQUE™ RIC 312 D WIDEX MOMENT™ , EVOKE™ PASSION WIDEX EVOKE™ , UNIQUE™	S	EASYWEAR Wired RIC Size: 0, 1, 2, 3, 4	INSTANT OPEN EAR-TIP XS, S, L TULIP EAR-TIP M, L DOUBLE EAR-TIP XS, S, M, L ROUND EAR-TIP, one-vent/two-vent; S, M, L CUSTOM OPEN EAR-TIP, TRANSPARENT/BEIGE HARD EAR-TIP, TRANSPARENT/BEIGE SOFT EAR-TIP, TRANSPARENT HARD/SOFT EXTENDED EAR-TIP, TRANSPARENT MODULAR/EMBEDDED EAR-TIP, TRANSPARENT
MODERATE	RIC 10 WIDEX MOMENT™ , EVOKE™ FUSION 2 WIDEX EVOKE™ FUSION WIDEX EVOKE™ , UNIQUE™ RIC 312 D WIDEX MOMENT™ , EVOKE™ PASSION WIDEX EVOKE™ , UNIQUE™	M	EASYWEAR Wired RIC Size: 0, 1, 2, 3, 4	INSTANT TULIP EAR-TIP M, L DOUBLE EAR-TIP XS, S, M, L ROUND EAR-TIP, one-vent/two-vent; S, M, L CUSTOM OPEN EAR-TIP, TRANSPARENT/BEIGE HARD EAR-TIP, TRANSPARENT/BEIGE SOFT EAR-TIP, TRANSPARENT HARD/SOFT EXTENDED EAR-TIP, TRANSPARENT MODULAR/EMBEDDED EAR-TIP, TRANSPARENT
SEVERE	RIC 312 D WIDEX MOMENT™ , EVOKE™ FUSION 2 WIDEX EVOKE™ FUSION WIDEX EVOKE™ , UNIQUE™	P	EASYWEAR Wired RIC Size: 0, 1, 2, 3, 4	INSTANT TULIP EAR-TIP M, L DOUBLE EAR-TIP XS, S, M, L ROUND EAR-TIP, one-vent/two-vent; S, M, L CUSTOM HARD EAR-TIP, TRANSPARENT/BEIGE SOFT EAR-TIP, TRANSPARENT HARD/SOFT EXTENDED EAR-TIP, TRANSPARENT MODULAR/EMBEDDED EAR-TIP, TRANSPARENT
SEVERE	RIC 312 D WIDEX MOMENT™ , EVOKE™ FUSION 2 WIDEX EVOKE™ FUSION WIDEX EVOKE™ , UNIQUE™	HP	EarWire RITE Size: -2, -1, 0, 1, 2, 3, 4, 5 Wired RITE Size: 0, 1, 2, 3, 4, 5	CUSTOM HARD EARMOULD*, TRANSPARENT/BEIGE MODULAR/EMBEDDED EAR-TIP, TRANSPARENT SOFT EARMOULD, TRANSPARENT

*Two options: Regular sound bore and Output Extender.

BTE MODELS

HEARING LOSS	MODELS/FAMILY	TUBE	SIZE	CUSTOM/INSTANT
MILD TO SEVERE	BTE 312 WIDEX MOMENT™ FASHION MINI WIDEX EVOKE™ , UNIQUE™ FASHION WIDEX EVOKE™ , UNIQUE™	EASYWEAR THIN TUBE 0.9 EASYWEAR THIN TUBE 1.4	EASYWEAR THIN TUBE Size: 0, 1, 2, 3, 4	INSTANT OPEN EAR-TIP XS, S, L TULIP EAR-TIP M, L DOUBLE EAR-TIP XS, S, M, L ROUND EAR-TIP, one-vent/two-vent; S, M, L CUSTOM OPEN EAR-TIP, TRANSPARENT/BEIGE HARD EAR-TIP, TRANSPARENT/BEIGE SOFT EAR-TIP, TRANSPARENT HARD/SOFT EXTENDED EAR-TIP, TRANSPARENT
MILD TO SEVERE / PROFOUND	FASHION MINI WIDEX EVOKE™ , UNIQUE™ FASHION WIDEX EVOKE™ , UNIQUE™ FASHION POWER WIDEX EVOKE™ , UNIQUE™ ME-SP* MENU™	CLASSIC MOULD	CLASSIC MOULD	CLASSIC MOULD OPEN EARMOULD HARD EARMOULD SOFT EARMOULD

* ME-SP is for PROFOUND hearing loss.

EASYWEAR INSTANT EAR-TIPS



CUSTOM EAR-TIPS



V.2 EASYWEAR WIRED RECEIVERS



V.2 EASYWEAR THIN TUBES



EASYWEAR INSTANT EAR-TIPS



CUSTOM EAR-TIPS



EASYWEAR WIRED RECEIVERS



EASYWEAR THIN TUBES



RITE WIRES



EARPIECE OVERVIEW

- EASYWEAR: FUSION 2, FUSION, PASSION, SUPER, FASHION, FASHION MINI AND FASHION POWER

WIDEX **MOMENT™** RIC 10, WIDEX **MOMENT™** RIC 312 D
WIDEX **MOMENT™** BTE 312

- V.2 EASYWEAR: WIDEX **MOMENT™** mRIC R D,
WIDEX **EVOKE™** BTE 13 D, WIDEX **MOMENT™** BTE 13 D

NEW!

UPDATED WITH
WIDEX **MOMENT™**
BTE 13 D
AND
BTE 312