






















































EARPIECE OVERVIEW

V.2 RIC MODELS

MODELS/FAMILY	HEARING LOSS	RECEIVER	SIZE	CUSTOM/INSTANT
 sRIC R D WIDEX MOMENT™	MILD TO MODERATE	 V.2 RIC M-Receiver	Size: 0, 1, 2 ,3, 4	INSTANT <ul style="list-style-type: none"> OPEN EAR-TIP <i>XS, S, L</i> OPEN ROUND EAR-TIP <i>S, M, L</i> ROUND EAR-TIP, one-vent <i>S, M, L</i> SLEEVE VENTED EAR-TIP <i>XS, S, M, L</i> SLEEVE POWER EAR-TIP <i>XS, S, M, L</i> CUSTOM <ul style="list-style-type: none"> OPEN EAR-TIP, <i>TRANSPARENT/BEIGE</i> SOFT EAR-TIP, <i>TRANSPARENT</i> HARD EAR-TIP, <i>TRANSPARENT</i> SOFT EXTENDED EAR-TIP, <i>TRANSPARENT</i> HARD EXTENDED EAR-TIP, <i>TRANSPARENT</i> EMBEDDED EAR-TIP, <i>TRANSPARENT</i>
				INSTANT <ul style="list-style-type: none"> ROUND EAR-TIP, one-vent <i>S, M, L</i> SLEEVE VENTED EAR-TIP <i>XS, S, M, L</i> SLEEVE POWER EAR-TIP <i>XS, S, M, L</i> CUSTOM <ul style="list-style-type: none"> OPEN EAR-TIP, <i>TRANSPARENT/BEIGE</i> SOFT EAR-TIP, <i>TRANSPARENT</i> HARD EAR-TIP, <i>TRANSPARENT</i> SOFT EXTENDED EAR-TIP, <i>TRANSPARENT</i> HARD EXTENDED EAR-TIP, <i>TRANSPARENT</i> EMBEDDED EAR-TIP, <i>TRANSPARENT</i>
 mRIC R D WIDEX MAGNIFY™	SEVERE	 V.2 RIC P-Receiver	Size: 0, 1, 2, 3, 4	INSTANT <ul style="list-style-type: none"> ROUND EAR-TIP, one-vent <i>S, M, L</i> SLEEVE VENTED EAR-TIP <i>XS, S, M, L</i> SLEEVE POWER EAR-TIP <i>XS, S, M, L</i> CUSTOM <ul style="list-style-type: none"> OPEN EAR-TIP, <i>TRANSPARENT/BEIGE</i> SOFT EAR-TIP, <i>TRANSPARENT</i> HARD EAR-TIP, <i>TRANSPARENT</i> SOFT EXTENDED EAR-TIP, <i>TRANSPARENT</i> HARD EXTENDED EAR-TIP, <i>TRANSPARENT</i> EMBEDDED EAR-TIP, <i>TRANSPARENT</i>
		 V.2 RIC HP-Receiver	Size: 0, 1, 2, 3, 4	CUSTOM <ul style="list-style-type: none"> EMBEDDED EAR-TIP, <i>TRANSPARENT</i>

V.2 THIN TUBES

MODELS/FAMILY	HEARING LOSS	TUBE	SIZE	CUSTOM/INSTANT
 BTE 312 WIDEX MOMENT™ , MAGNIFY™	MILD TO SEVERE	 V.2 THIN TUBE 0.9	Size: 0, 1, 2 ,3, 4	INSTANT <ul style="list-style-type: none"> OPEN EAR-TIP <i>XS, S, L</i> OPEN ROUND EAR-TIP <i>S, M, L</i> ROUND EAR-TIP, one-vent <i>S, M, L</i> SLEEVE VENTED EAR-TIP <i>XS, S, M, L</i> SLEEVE POWER EAR-TIP <i>XS, S, M, L</i> CUSTOM <ul style="list-style-type: none"> OPEN EAR-TIP, <i>TRANSPARENT/BEIGE</i> SOFT EAR-TIP, <i>TRANSPARENT</i> HARD EAR-TIP, <i>TRANSPARENT</i> SOFT EXTENDED EAR-TIP, <i>TRANSPARENT</i> HARD EXTENDED EAR-TIP, <i>TRANSPARENT</i>
		 V.2 THIN TUBE 1.4		CUSTOM <ul style="list-style-type: none"> OPEN EAR-TIP, <i>TRANSPARENT/BEIGE</i> SOFT EAR-TIP, <i>TRANSPARENT</i> HARD EAR-TIP, <i>TRANSPARENT</i> SOFT EXTENDED EAR-TIP, <i>TRANSPARENT</i> HARD EXTENDED EAR-TIP, <i>TRANSPARENT</i>



































EARPIECE OVERVIEW

RIC/RITE MODELS

MODELS/FAMILY	HEARING LOSS	RECEIVER	SIZE	CUSTOM/INSTANT
 RIC 10 WIDEX MOMENT™ , MAGNIFY™	MILD	 Wired S-RIC		INSTANT <ul style="list-style-type: none"> OPEN EAR-TIP <i>XS, S, L</i> OPEN ROUND EAR-TIP <i>S, M, L</i> ROUND EAR-TIP, one-vent <i>S, M, L</i> SLEEVE VENTED EAR-TIP <i>XS, S, M, L</i> SLEEVE POWER EAR-TIP <i>XS, S, M, L</i> CUSTOM <ul style="list-style-type: none"> OPEN EAR-TIP, <i>TRANSPARENT/BEIGE</i> SOFT EAR-TIP, <i>TRANSPARENT</i> HARD EAR-TIP, <i>TRANSPARENT</i> SOFT EXTENDED EAR-TIP, <i>TRANSPARENT</i> HARD EXTENDED EAR-TIP, <i>TRANSPARENT</i> EMBEDDED EAR-TIP, <i>TRANSPARENT</i>
				Size: 0, 1, 2 ,3, 4
 PASSION WIDEX EVOKE™	MODERATE	 Wired M-RIC		INSTANT <ul style="list-style-type: none"> OPEN EAR-TIP, <i>TRANSPARENT/BEIGE</i> SOFT EAR-TIP, <i>TRANSPARENT</i> HARD EAR-TIP, <i>TRANSPARENT</i> SOFT EXTENDED EAR-TIP, <i>TRANSPARENT</i> HARD EXTENDED EAR-TIP, <i>TRANSPARENT</i> EMBEDDED EAR-TIP, <i>TRANSPARENT</i> CUSTOM <ul style="list-style-type: none"> OPEN EAR-TIP, <i>TRANSPARENT/BEIGE</i> SOFT EAR-TIP, <i>TRANSPARENT</i> HARD EAR-TIP, <i>TRANSPARENT</i> SOFT EXTENDED EAR-TIP, <i>TRANSPARENT</i> HARD EXTENDED EAR-TIP, <i>TRANSPARENT</i> EMBEDDED EAR-TIP, <i>TRANSPARENT</i>
				Size: 0, 1, 2, 3, 4
 RIC 312 D WIDEX MOMENT™ , EVOKE™ , MAGNIFY™	SEVERE	 Wired P-RIC		INSTANT <ul style="list-style-type: none"> ROUND EAR-TIP, one-vent <i>S, M, L</i> SLEEVE VENTED EAR-TIP <i>XS, S, M, L</i> SLEEVE POWER EAR-TIP <i>XS, S, M, L</i> CUSTOM <ul style="list-style-type: none"> OPEN EAR-TIP, <i>TRANSPARENT/BEIGE</i> SOFT EAR-TIP, <i>TRANSPARENT</i> HARD EAR-TIP, <i>TRANSPARENT</i> SOFT EXTENDED EAR-TIP, <i>TRANSPARENT</i> HARD EXTENDED EAR-TIP, <i>TRANSPARENT</i>
				Size: 0, 1, 2, 3, 4
 FUSION 2 WIDEX EVOKE™	SEVERE	 HP		CUSTOM <ul style="list-style-type: none"> HARD EARMOULD*, <i>TRANSPARENT/BEIGE</i> MODULAR/EMBEDDED EAR-TIP, <i>TRANSPARENT</i> SOFT EARMOULD, <i>TRANSPARENT</i>
				EarWire RITE Size: -2, -1, 0, 1, 2, 3, 4, 5
 RIC 312 D WIDEX MOMENT™ , MAGNIFY™	SEVERE			Wired RITE Size: 0, 1, 2, 3, 4, 5

*Two options: Regular sound bore and Output Extender.

BTE MODELS

MODELS/FAMILY	HEARING LOSS	TUBE	SIZE	CUSTOM/INSTANT
 FASHION MINI WIDEX EVOKE™	MILD TO SEVERE	 EASYWEAR THIN TUBE 0.9		INSTANT <ul style="list-style-type: none"> OPEN EAR-TIP <i>XS, S, L</i> OPEN ROUND EAR-TIP <i>S, M, L</i> ROUND EAR-TIP, one-vent <i>S, M, L</i> SLEEVE VENTED EAR-TIP <i>XS, S, M, L</i> SLEEVE POWER EAR-TIP <i>XS, S, M, L</i> CUSTOM <ul style="list-style-type: none"> OPEN EAR-TIP, <i>TRANSPARENT/BEIGE</i> SOFT EAR-TIP, <i>TRANSPARENT</i> HARD EAR-TIP, <i>TRANSPARENT</i> SOFT EXTENDED EAR-TIP, <i>TRANSPARENT</i> HARD EXTENDED EAR-TIP, <i>TRANSPARENT</i>
				Size: 0, 1, 2 ,3, 4
 FASHION WIDEX EVOKE™		 EASYWEAR THIN TUBE 1.4		INSTANT <ul style="list-style-type: none"> ROUND EAR-TIP, one-vent <i>S, M, L</i> SLEEVE VENTED EAR-TIP <i>XS, S, M, L</i> SLEEVE POWER EAR-TIP <i>XS, S, M, L</i> CUSTOM <ul style="list-style-type: none"> OPEN EAR-TIP, <i>TRANSPARENT/BEIGE</i> SOFT EAR-TIP, <i>TRANSPARENT</i> HARD EAR-TIP, <i>TRANSPARENT</i> SOFT EXTENDED EAR-TIP, <i>TRANSPARENT</i> HARD EXTENDED EAR-TIP, <i>TRANSPARENT</i>
				Size: 0, 1, 2, 3, 4
 BTE 312 WIDEX MOMENT™ , MAGNIFY™	MILD TO SEVERE / PROFOUND	 BTE 13 D WIDEX MOMENT™ , MAGNIFY™		INSTANT <ul style="list-style-type: none"> OPEN EARMOULD CUSTOM <ul style="list-style-type: none"> HARD EARMOULD
				INSTANT <ul style="list-style-type: none"> OPEN EARMOULD CUSTOM <ul style="list-style-type: none"> HARD EARMOULD
 FASHION MINI WIDEX EVOKE™		 FASHION WIDEX EVOKE™		INSTANT <ul style="list-style-type: none"> OPEN EARMOULD CUSTOM <ul style="list-style-type: none"> HARD EARMOULD
				INSTANT <ul style="list-style-type: none"> OPEN EARMOULD CUSTOM <ul style="list-style-type: none"> HARD EARMOULD

* ME-SP is for PROFOUND hearing loss.

EASYWEAR - RECOMMENDED INSTANT EAR-TIPS

V.2 RIC M AND P-RECEIVERS

V.2 THIN TUBE 0.9 AND 1.4

OPEN - SIZE: XS, S, L

OPEN ROUND - SIZE: S, M, L

ROUND - SIZE: S, M, L

ONE-VENT

SLEEVE VENTED: XS, S, M, L

SLEEVE POWER: XS, S, M, L

CUSTOM EAR-TIPS

V.2 RIC M-RECEIVERS	V.2 RIC M AND P-RECEIVERS					M, P
V.2 THIN TUBE 0.9 AND 1.4						
OPEN EAR-TIP	SOFT EAR-TIP	HARD EAR-TIP	SOFT EXTENDED EAR-TIP	HARD EXTENDED EAR-TIP	EMBEDDED EAR-TIP	

TRANSPARENT/BEIGE

V.2 EASYWEAR WIRED RECEIVERS

V.2 RIC M AND P-RECEIVERS

SIZE: 0, 1, 2, 3, 4
LEFT/RIGHT

V.2 RIC HP RECEIVER

SIZE: 0, 1, 2, 3, 4
LEFT/RIGHT
EMBEDDED TRANSPARENT ONLY

V.2 EASYWEAR THIN TUBES

V.2 THIN TUBE 0.9 AND 1.4

SIZE: 0, 1, 2, 3, 4
LEFT/RIGHT

EASYWEAR - RECOMMENDED INSTANT EAR-TIPS

WIRED S, M, P-RECEIVERS

EASYWEAR THIN TUBE 0.9 AND 1.4

OPEN - SIZE: XS, S, L

OPEN ROUND - SIZE: S, M, L

ROUND - SIZE: S, M, L

ONE-VENT

SLEEVE VENTED: XS, S, M, L

SLEEVE POWER: XS, S, M, L

CUSTOM EAR-TIPS

WIRED S, M	WIRED S, M, P					M, P
EASYWEAR THIN TUBE 0.9 AND 1.4						
OPEN EAR-TIP	SOFT EAR-TIP	HARD EAR-TIP	SOFT EXTENDED EAR-TIP	HARD EXTENDED EAR-TIP	EMBEDDED EAR-TIP	

TRANSPARENT/BEIGE

EASYWEAR WIRED RECEIVERS

EASYWEAR WIRED S, M AND P-RECEIVERS

SIZE: 0, 1, 2, 3, 4
LEFT/RIGHT

EASYWEAR THIN TUBES

EASYWEAR THIN TUBE 0.9 AND 1.4

SIZE: 0, 1, 2, 3, 4
LEFT/RIGHT

RITE WIRES

RITE WIRES	CUSTOM EARMOULDS	
SIZE: -2, -1, 0, 1, 2, 3, 4, 5	RITE HP, SP HARD EARMOULD	WIRED RITE HP, SP SOFT EARMOULD

EARPIECE OVERVIEW

V.2 EASYWEAR
EASYWEAR



UPDATED
2023 VERSION

